

arlington

clinical recliner


rt3281

dimensions
36w . 39d . 46h



healthcare furnishings

|durfold|



arlington clinical recliner

rt3281

3" diameter, locking, non-marring front and rear casters for ease of steering and mobility

infinite position reclining system to accommodate all users and circumstances including the fowler, supine and trendelenberg positions

heavy duty steel handle for footpiece and trendelenberg position activation to endure consistent use

able to be reclined to a horizontal position for 72" long sleeping surface

optional recessed sliding foot tray for ease of transporting patients

additional options (for all of our recliners) include push bar, tray table, iv pole and attachment, foley bag holder, user restraint attachments, patient alarm storage pocket

certified bariatric testing results"

GSA Cyclic Increasing Back Load Test

Passed Max Load = 175 lbs.

Passed Max Cycles = 125,000

GSA Cyclic Vertical Load Test on One Arm

Passed Max Load = 250 lbs.

Passed Max Cycles = 100,000

GSA Cyclic Side Thrust Load Test on Arms

Passed Max Load = 200 lbs.

Passed Max Cycles = 175,000

GSA Cyclic Vertical Load Test on Seats, Bases and Casters

Passed Max Load = 1300 lbs.

Passed Max Cycles = 300,000

GSA Test on Front Stability

Passed Max Load = 400 lbs.

GSA Test on Back Stability

Passed Max Load = 2002.2 in-lbs.

* durfold reserves the right to make changes to products or information at any time without notice

	recliner
w	36*
d	39
h	46
sw	23
sd	21
sh	22
ah	28
sl	72

approximate dimensions given

* includes activator levers. overall width of arms = 30"

durfold corporation
 102 upton drive, jackson, ms 39209
 p.800.345.6849 f.601.922.6244
 durfold.com

|durfold|